



## **STALHAM HOISTS**

Lifting & Sit-stand

Versatile, portable hoists for patient positioning

# STALHAM 150 LIFTING HOIST



## Structure

Epoxy coated steel frame.

Complete lifting hoist can be assembled in just a few minutes.

Utilises electric lifting arm.



## Comfort

A range of comfortable slings are available to fit patient requirements.

Ergonomic handles.



## Ease Of Handling

Wired remote control for user controls in any position around the lift.

Capable of lowering patient to the ground.

The base frame has a very low construction to allow for very good under-ride possibilities.

Undercarriage is expandable as needed by means of foot pedal.

Mounted on four swivelling castors for portability. Each rear castor is braked.



## Safety

Features easy to access emergency stop button.

Lift can allow for emergency mechanical lowering if required.

Additional battery backup.

Medical grade electric lifting arm.

# STALHAM 150 SIT-STAND LIFT



## Structure

Epoxy coated steel frame.

Complete lifting hoist can be assembled in just a few minutes.

Utilises electric lifting arm.



## Comfort

Leg straps ensure comfort and stability for patient when lifting.

180° grab handles support patient.

Angled foot rests secure feet in place when in use.



## Ease Of Handling

Wired remote control for user controls in any position around the lift.

The base frame has a very low construction to allow for very good under-ride possibilities.

Undercarriage is expandable as needed by means of foot pedal.

Mounted on four swivelling castors for portability. Each rear castor is braked.



## Safety

Features easy to access emergency stop button.

Lift can allow for emergency mechanical lowering if required.

Additional battery backup.

Medical grade electric lifting arm.

## Specifications

# STALHAM 150 LIFTING HOIST

### Hoist: 150 LIFTING HOIST

Lift height (min.)	55 cm
Lift height (max.)	174 cm
Underpin height	14 cm
Device weight	43 kg
Max. user weight	150 kg

### Accessories:

Toilet sling (optional velcro fastener version)  
Approved up to 250 kg weight



Universal sling (S, M/L, XL/XXL)  
Approved up to 250 kg weight



Universal sling with headrest (M/L)  
Approved up to 300 kg weight



Bariatric sling  
Approved up to 300 kg weight



Net/bath sling (S, M/L, XL/XXL)  
Approved up to 250 kg weight




Wide Hanger  
4-point spreader bar 56 x 33 x 14cm



Horizontal sling (optional split leg)  
Approved up to 250 kg weight



Hoist Scale PW500  
Approved to class III, MPG class I



## Specifications

# STALHAM 150 SIT-STAND HOIST

### Hoist: 150 SIT-STAND LIFT

Lift height (min.)	84 cm
Lift height (max.)	150 cm
Underpin height	14 cm
Device weight	53 kg
Max. user weight	150 kg

### Accessories:

#### Back support sling

Approved up to 250 kg weight



#### Calf support belt



[www.medi-plinth.co.uk](http://www.medi-plinth.co.uk)

7-11 Holywells Road, Ipswich, Suffolk IP3 0DL – Call: +44 (0) 1473 212010