



STALHAM BATHE TROLLEYS

Height adjustable shower trolleys

Versatile, comfortable hygiene solutions

Hydraulic Shower Trolley

STALHAM BATHE 20



Structure

The chassis and sub-frames are made of powder coated stainless steel.

Four different models of the Bathe 20 are available, to ensure a suitable size for most care environments.

Attractive design, with mattress available in a selection of colours.

Lift system supports weights of up to 200 kg.



Comfort

Padded mattress included for patient comfort, is easily removable for cleaning.

Fasteners/straps are available as accessories, for attachment on mattresses to siderails.

Trolley drains at both ends allowing the user person to be positioned in either direction.

The stretcher can easily be tilted in both directions to the Trendelenburg position.

Hydraulic height adjustment between 60–90cm



Ease Of Handling

The stretcher has rubber bumper wheels on all four corners to protect walls and doorways.

Floor clearance at 160 mm to enable floor lift transfers.

Foldable side rails for efficient and safe patient transfer and transport. The side rails can be pressed together when travelling through narrow doorways.

Dual drainage offers easy change of designated head end.



Controls

Foot operated hydraulic pedals allow for control of height adjustment.

The stretcher lowers smoothly when pedal is pressed and stops immediately when pedal is released.

Manual tilt adjustment can be accessed from both sides of the trolley.

Electric Shower Trolley

STALHAM BATHE 30



Structure

The chassis and sub-frames are made of powder coated stainless steel.

Three lengths are available, to ensure a suitable size for most care environments.

Attractive design, with mattress available in a selection of colours.

Lift system supports weights of up to 200 kg.



Comfort

Padded mattress included for patient comfort, is easily removable for cleaning.

Fasteners/straps are available as accessories, for attachment on mattresses to siderails.

Trolley drains at both ends allowing the user person to be positioned in either direction.

The stretcher can easily be tilted in both directions to the Trendelenburg position.

Electric height adjustment between 60–90cm



Ease Of Handling

Includes an on-board battery solution as standard. Also available with a removable battery option.

The stretcher has rubber bumper wheels on all four corners to protect walls and doorways.

Floor clearance at 160 mm to enable floor lift transfers.

Foldable side rails for efficient and safe patient transfer and transport. The side rails can be pressed together when travelling through narrow doorways.

Dual drainage offers easy change of designated head end.



Controls

Audible tone when battery is running low.

Foot controller as standard, allows for control of electric height adjustment.

Emergency stop button located on battery.

Manual tilt adjustment can be accessed from both sides of the trolley.

Specifications

STALHAM BATHE 20 & 30

Hydraulic Trolley: BATHE 20

Trolley Size	160 x 66cm
Patient area	146 x 56cm
Trolley Weight	88kg

Electric Trolley: BATHE 30

Trolley length	160 x 66cm
Patient area	146 x 56cm
Trolley Weight	88kg

Short Hydraulic Trolley: BATHE 20S

Trolley length	190 x 66cm
Patient area	176 x 56cm
Trolley Weight	80kg

Short Electric Trolley: BATHE 30S

Trolley length	190 x 66cm
Patient area	176 x 56cm
Trolley Weight	80kg

Long Hydraulic Trolley: BATHE 20L

Trolley length	210 x 66cm
Patient area	196 x 56cm
Trolley Weight	92kg

Long Electric Trolley: BATHE 30L

Trolley length	210 x 66cm
Patient area	196 x 56cm
Trolley Weight	92kg

Wide Hydraulic Trolley: BATHE 20W

Trolley length	195 x 70cm
Patient area	178 x 64cm
Trolley Weight	91kg

Electric features:

Charger input	110-240V / 50-60Hz
Supply voltage	24 V DC



www.medi-plinth.co.uk

7-11 Holywells Road, Ipswich, Suffolk IP3 0DL – Call: +44 (0) 1473 212010