

## Ward Furniture

# IV POLES / DRIP STANDS



### Features

These easy-assembly drip stands are manufactured from medical grade stainless steel and high-quality injection moulded parts. All IV poles come with a quick clip drop catcher. The weighted base models include a sturdy handle for help with moveability.

The five-legged bases with swivel castors ensure excellent stability and manoeuvrability. The adjustable height range with simple locking mechanism are perfect for most care home or hospital environments.

- To reduce transport costs and unnecessary packaging IV poles are flat packed.
- Drop catcher with quick clip.
- Two hook or four hook.
- Handle with included with weighted base models.
- Two hook wall mounted version also available.
- Five ESD (Electrostatic dissipative) swivel castors, two of which are locking.
- Drips Stands comply to BS ISO 15375:2010

### IV Poles: MOBILE / WALL MOUNTED

Flat pack 2-hook mobile: **IVDST001**

Flat pack 2-hook mobile, weighted base: **IVDST002**

Flat pack 4-hook mobile: **IVDST003**

Flat pack 4-hook mobile, weighted base: **IVDST004**

2-hook wall mounted: **IVDST005**

|                   |            |
|-------------------|------------|
| Base diameter     | Ø 71 cm    |
| Adjustable height | 129–204 cm |